

RoyIntelli Technology Co., Ltd.
RI-CHN2200 | CHN Elemental Analyzer
Category: Proximate & Ultimate Analysis
Buyer datasheet | Version 1.0 | 2026-09-24

Purpose and sample types

- Coal, coke, biomass, SRF, petroleum and chemicals
- ASTM D5373-14; GB/T 30733-2014
- C/H/N analysis in about 300 s per sample
- C 0.005-100%; H 0.05-25%; N 0.01-50%; 36-position autosampler
- Infrared detection for C/H; thermal conductivity detection for N
- Factory documents, English manual and remote commissioning are available according to the confirmed order

Ordering checklist

- Confirm sample matrix, daily throughput, reporting standard, voltage and plug format.
- Request the model-specific quotation, acceptance protocol and available evidence pack.
- Final specifications, documents, lead time, warranty and payment terms are confirmed in the Proforma Invoice.

Documents and support

- English operating documentation and remote commissioning are available according to the confirmed order.
- CRM reports, calibration records, conformity documents and reference material are supplied or made available only when applicable and verifiable for the ordered model.

Contact

service@royintelli.com | +86 17808359107 | WhatsApp +852 57166955
www.royintelli.com | Changsha, Hunan, China

This datasheet is a procurement summary, not a certificate. Confirm method compatibility and all commercial details before purchase.