

RoyIntelli Technology Co., Ltd.
RI-MAC67 | Automatic Proximate Analyzer (TGA)
Category: Proximate & Ultimate Analysis
Buyer datasheet | Version 1.0 | 2026-09-24

Purpose and sample types

- Coal, coke, biomass, solid waste
- ASTM D7582; GB/T 212-2008; GB/T 211-2007; GB/T 30729-2014; YS/T 587; HJ 1024-2019
- 20 samples; 90-120 min for M, A and V
- Ambient-1050 C; temperature control +/-1 C; weighing resolution 0.0001 g
- Dual-furnace operation; built-in thermobalance; silent air pump
- Factory documents, English manual and remote commissioning are available according to the confirmed order

Ordering checklist

- Confirm sample matrix, daily throughput, reporting standard, voltage and plug format.
- Request the model-specific quotation, acceptance protocol and available evidence pack.
- Final specifications, documents, lead time, warranty and payment terms are confirmed in the Proforma Invoice.

Documents and support

- English operating documentation and remote commissioning are available according to the confirmed order.
- CRM reports, calibration records, conformity documents and reference material are supplied or made available only when applicable and verifiable for the ordered model.

Contact

service@royintelli.com | +86 17808359107 | WhatsApp +852 57166955
www.royintelli.com | Changsha, Hunan, China

This datasheet is a procurement summary, not a certificate. Confirm method compatibility and all commercial details before purchase.