

RoyIntelli Technology Co., Ltd.
RI-MW2200 | Automatic Moisture Analyzer
Category: Proximate & Ultimate Analysis
Buyer datasheet | Version 1.0 | 2026-09-24

Purpose and sample types

- Coal and coke
- GB/T 211-2007; GB/T 212-2008; GB/T 483-2007
- 22 samples; classic 25-45 min or fast 15-25 min
- Light-wave heating; optional nitrogen purge; power <2.5 kW
- Automatic weighing; heat dissipation; power-loss data memory
- Factory documents, English manual and remote commissioning are available according to the confirmed order

Ordering checklist

- Confirm sample matrix, daily throughput, reporting standard, voltage and plug format.
- Request the model-specific quotation, acceptance protocol and available evidence pack.
- Final specifications, documents, lead time, warranty and payment terms are confirmed in the Proforma Invoice.

Documents and support

- English operating documentation and remote commissioning are available according to the confirmed order.
- CRM reports, calibration records, conformity documents and reference material are supplied or made available only when applicable and verifiable for the ordered model.

Contact

service@royintelli.com | +86 17808359107 | WhatsApp +852 57166955
www.royintelli.com | Changsha, Hunan, China

This datasheet is a procurement summary, not a certificate. Confirm method compatibility and all commercial details before purchase.